**WARNING:** TO PREVENT SERIOUS OR FATAL INJURY: FOR USE ONLY BY PERSONS WEIGHING LESS THAN 300 POUNDS. Consult healthcare professional for advice and instructions specific to your needs before use. Do not jump or stand on seat. Contact CSI immediately at 1-800-882-0116 if any components are missing or damaged. Assemble in accordance with CSI's instructions. Do not make substitutions or modify this product in any way. Wipe seat dry after use. Check seat periodically after assembly to make sure all screws are tight. FOR QUESTIONS OR ASSISTANCE WITH INSTALLATION CONTACT CSI AT 1-800-882-0116.

**LIMITED LIFETIME WARRANTY:** If this product fails due to a defect in materials or workmanship at any time during the life of the product, Creative Specialties International, a division of Moen Incorporated (CSI) will replace it free of charge, postage-paid. Simply contact CSI at 1-800-882-0116 for replacement information. This warranty does not cover products which have been abused, altered, damaged, misused, cut or worn. TO THE EXTENT PERMITTED BY LAW, CSI DISCLAIMS ALL OTHER IMPLIED OR EXPRESS WARRANTIES INCLUDING ALL WARRANTIES OF MERCHANTABILITY AND/OR FITNESS FOR A PARTICULAR PURPOSE.

**LIMITATION OF LIABILITY:** THIS WARRANTY IS YOUR SOLE REMEDY AND CSI SHALL NOT BE LIABLE FOR ANY DAMAGES, WHETHER DIRECT, INDIRECT, INCIDENTAL, SPECIAL, CONSEQUENTIAL, EXEMPLARY, OR OTHERWISE, INCLUDING LOST REVENUES AND LOST PROFITS, ARISING OUT OF ANY THEORY OR RECOVERY, INCLUDING STATUTORY, CONTRACT OR TORT. NOTWITHSTANDING THE TERM OF ANY LIMITED OR IMPLIED WARRANTY, OR IN THE EVENT THAT ANY LIMITED WARRANTY FAILS OF ITS ESSENTIAL PURPOSE, IN NO EVENT WILL CSI'S ENTIRE LIABILITY EXCEED THE PURCHASE PRICE OF THIS PRODUCT. SOME STATES AND PROVINCES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES SO THE ABOVE LIMITATIONS OR EXCLUSIONS MAY NOT APPLY TO YOU. THIS LIMITED WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE AND PROVINCE TO PROVINCE.

For additional warranty information, visit [www.homecare.moen.com](http://www.homecare.moen.com).

The Home Care Adjustable Transfer Bench meets or exceeds the applicable requirements of the following standards / guidelines:

 The Americans with Disabilities Act (ADA)

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
INNOVATIVE DESIGNS  
for independent living™

## DN8065 ADJUSTABLE TRANSFER BENCH INSTALLATION INSTRUCTIONS

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Home Care is a registered trademark of Home Care Tim-BR-MART, licensed to Moen Incorporated. Home Care and Creative Specialties are registered trademarks of Moen Incorporated. Creative Specialties is a division of Moen Incorporated. © 2006

Home Care es una marca registrada de Home Care Tim-BR-MART, con licencia a Moen Incorporated. Home Care y Creative Specialties son marcas registradas de Moen Incorporated. Creative Specialties es una división de Moen Incorporated. © 2006

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
INNOVATIVE DESIGNS  
for independent living™

**DN8065**  
**ADJUSTABLE TRANSFER BENCH**  
**INSTALLATION INSTRUCTIONS**  
**SILLA DE TRANSFERENCIA AJUSTABLE**  
**INSTRUCCIONES DE INSTALACIÓN**  
**BANC DE TRANSFERT AJUSTABLE**  
**DIRECTIVES D'INSTALLATION**

CONTENTS:  
CONTENIDO:  
CONTENU:

		1 seat back 1 respaldo 1 dos de banc		Hardware: 11 1-1/2" bolts, 11 washers, 11 lock washers, 4 lock nuts	  
4 adjustable legs 4 patas ajustables 4 pattes ajustables		1 transfer bench seat bottom 1 fondo de asiento de silla de transferencia 1 banc de transfert (extérieur)		Ferretería: 11 pernos de 1-1/2" (3.81 cm), 11 arandelas, 11 arandelas de seguridad, 4 tuercas de seguridad	
2 support rails 2 barandas de apoyo 2 barres de support		1 large seat bottom 1 asiento grande 1 grand banc de transfert (intérieur)		Quincaillerie : 11 boulons de 1-1/2 po, 11 rondelles, 11 rondelles de blocage, 4 contre-écrous	

**WARNING! - FOLLOW THESE INSTALLATION INSTRUCTIONS CAREFULLY. PROPER INSTALLATION IS THE INSTALLER'S RESPONSIBILITY. FAILURE TO FOLLOW CORRECT INSTALLATION PROCEDURES CAN CAUSE THE SEAT TO BECOME LOOSE, WHICH CAN RESULT IN SERIOUS INJURY.**

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**Specifications:**

Fits inside most standard size tubs.  
Wipe seat dry after use. Check seat periodically after assembly to make sure all screws are tight.

**Dimensions:**

33-1/4" wide x 19-1/4" deep x 17" high (at lowest leg setting).

**Tools Needed:**

Standard size Phillips head screwdriver.

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
Determine direction chair must point to face water supply. Install seat bottom with large holes to the back of the transfer bench seat bottom. (figure 4)

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