


Regal 190 S.S. — Rev. 1 (12/01) Copyright © 2001 SLOAN VALVE COMPANY Printed in the U.S.A.

# Regal<sup>®</sup> Model Flushometer



## Description

Concealed Urinal Flushometer, for 11/4" back spud urinals.

### Flush Cycle

□ Model 190 Conventional (3.5 gpf/13.2 Lpf) □ Model 190-1.5 Water Saver (1.5 gpf/5.7 Lpf) □ Model 190-1.0 Low Consumption (1.0 gpf/3.8 Lpf)

### Specifications

Quiet, Concealed, Diaphragm Type, Rough Brass Urinal Flushometer for either left or right hand supply with the following features:

- ADA Compliant Metal Oscillating Non-Hold-Open Handle
- 1" I.P.S. Wheel Handle Bak-Chek™ Angle Stop
- Adjustable Tailpiece
- Vacuum Breaker
- Elbow Flush Connection and Spud Coupling for 1<sup>1</sup>/<sub>4</sub>" Concealed Back Spud
- Exposed Parts Chrome Plated
- High Copper, Low Zinc Brass Castings for Dezincification Resistance
- Non-Hold-Open Handle and No External Volume Adjustment to Ensure Water Conservation
- Low Consumption flush accuracy controlled by Para-Flo<sup>™</sup> Technology
- Handle Packing, Stop Seat and Vacuum Breaker molded from PERMEX<sup>™</sup> Rubber Compound for Chloramine Resistance

Valve Body, Cover, Tailpiece and Control Stop shall be in conformance with ASTM Alloy Classification for Semi-Red Brass. Valve shall be in compliance to the applicable sections of ASSE 1037, ANSI/ASME 112.19.6, and Military Specification V-29193.

### L Dimension

Specify the "L" Dimension for the proper length of the Handle Assembly and Flush Connection. The "L" Dimension is equal to the Wall Thickness (to the nearest whole inch) plus 23/4" (70 mm).

### Variations

- □ T 1½" Flush Connection
- WB
   For installation with Sloan Easy Access® Wall Box

   YB
   Sweat Solder Adapter

See Accessories Section of the Sloan catalog for details on these and other Flushometer variations.





This space for Architect/Engineer approval

The information contained in this document is subject to change without notice.



 Made in the U.S.A.

 SLOAN VALVE COMPANY • 10500 SEYMOUR AVE. • FRANKLIN PARK, IL. 60131

 Ph: 1-800-9-VALVE-9 or 1-847-671-4300 • Fax: 1-800-447-8329 or 1-847-671-4380

 http://www.sloanvalve.com